




# DynaClose®

with **ELASTIC** traction

**Pre-op Skin Expansion**  
**Topical Skin Closure**

Soft fabric stretches to prevent blistering and dries quickly after showering.

Silicone elastomer provides elastic traction and won't stick to tissue.



[www.dynaclose.com](http://www.dynaclose.com)

For sales or distribution in your area, contact:  
[sales@canica.com](mailto:sales@canica.com)

1D16347A

**CANICA®** USA/Can: 1.800.705.8312  
Global: +1.613.256.0350  
[www.canica.com](http://www.canica.com) [www.dynaclose.com](http://www.dynaclose.com)

PRINTED IN CANADA


© CANICA DESIGN INC.

# DynaClose®

**Pre-op Skin Expansion**



Before defect excision



- ▶ Topically applied in office 1-2 weeks pre-op.
- ▶ Patient managed until day of surgery.
- ▶ Creates a skin depot with no surgical intervention.
- ▶ Enables primary closure after excision surgery.
- ▶ Facilitates a low-tension environment for faster healing with less pain and reduced scarring<sup>1</sup>.

After




1. *Healing Time of Radial Forearm Free Flap Donor Sites After Pre-operative Tissue Expansion: A Randomized Controlled Trial.* Bonaparte JP, Corsten M, Allen M. Ottawa Hospital, Ottawa, Canada. Presented May, 2010

# DynaClose®

**Topical Skin Closure**



Before



- ▶ Achieves delayed primary closure without sutures, staples, grafts or negative pressure.
- ▶ Easy and painless application is ideal for children and seniors.
- ▶ Maintains primary closure in a low-tension healing environment.
- ▶ Controls infected wound edges, then closes them when ready.

After



2. *Saving the failing wound.* Bell MSG. Canadian Journal of Plastic Surgery. 2009 Fall; 17(3):e15-e16

**CANICA®**